


# Hydrofera Blue READY and Hydrofera Blue READY-Transfer Dressings




## Bioburden and Exudate Control at Your Fingertips

- Available with or without film backing
- Antibacterial properties to help manage bioburden
  - No hydration required prior to application
  - Up to 7-day wear time

*Hydrofera*<sup>®</sup>

# When the Management of Bioburden Is a Priority



## Broad-spectrum antibacterial protection<sup>1</sup>

- Provided by two organic pigments, methylene blue and gentian violet
- Noncytotoxic<sup>2</sup>

## Absorption

- Absorbs and retains bacteria-laden exudate

Available with or without waterproof film backing

## Hydrofera Blue READY Dressings Provide Broad-Spectrum Antibacterial Protection<sup>1</sup>

### Bacteria

- *Methicillin-resistant Staphylococcus aureus (MRSA)*
- *Staphylococcus aureus*
- *Staphylococcus epidermidis (coagulase negative)*
- *Pseudomonas aeruginosa*
- *Pseudomonas fluorescens*
- *Escherichia coli*
- *Enterococcus faecalis*
- *Vancomycin-resistant Enterococcus (VRE)*
- *Streptococcus pneumoniae*
- *Streptococcus pyogenes*
- *Klebsiella pneumoniae*
- *Proteus mirabilis*
- *Proteus vulgaris*
- *Enterobacter aerogenes*
- *Enterobacter cloacae*
- *Bacillus subtilis*
- *Yersinia enterocolitica*
- *Serratia marcescens*

### Yeast

- *Candida albicans*
- *Candida krusei*
- *Candida glabrata*



## Hydrofera Blue READY



- Waterproof film backing
- No secondary dressing required
- Exudate is absorbed and retained
- Up to 7-day wear time

## Hydrofera Blue READY-Transfer




- No film backing increases exudate management options
- Exudate is transferred into secondary absorbent dressing
- Suitable for use on partial and full thickness wounds
- Up to 7-day wear time

### Performance Data

#### Absorption Comparison

**Hydrofera Blue READY** foam dressing has absorption comparable to other antimicrobial foam dressings, based on in vitro testing<sup>3</sup>




(Error bars represent one standard deviation)

### Performance Data

#### In Vitro Fluid Transfer Study

**Hydrofera Blue READY-Transfer** has a high rate of exudate transfer to a secondary dressing<sup>4</sup>



(Error bars represent one standard deviation)

## Hydrofera Blue **READY-Transfer** antibacterial dressings

Stock No	Product Size	Quantity/Box	HCPCS*
<b>Standard Dressings</b>			
HBRT2525	2.5" x 2.5" (6.4 cm x 6.4 cm)	10	A6209
HBRT4050	4" x 5" (10.2 cm x 12.7 cm)	10	A6210
HBRT8080	8" x 8" (20.3 cm x 20.3 cm)	10	A6211



No hydration required

## Hydrofera Blue **READY** antibacterial dressings


Stock No	Product Size	Quantity/Box	HCPCS*
<b>Standard Dressings</b>			
HBRS2520	2.5" x 2.5" (6.4 cm x 6.4 cm)	10	A6209
HBRS4520	4" x 5" (10.2 cm x 12.7 cm)	10	A6210
HBRS8820	8" x 8" (20.3 cm x 20.3 cm)	10	A6211



No hydration required

## Hydrofera Blue **CLASSIC** antibacterial dressings

Stock No	Product Size	Quantity/Box	HCPCS*
<b>Standard Dressings</b>			
HB2214	2" x 2" (5 cm x 5 cm)	10	A6209
HB4414	4" x 4" (10 cm x 10 cm)	10	A6209
HB6614	6" x 6" (15 cm x 15 cm)	10	A6210
HBT0906	9-mm Tunneling Dressing (1.2 g)	10 pair	A6215
<b>Heavy Drainage Dressings</b>			
HBHD4450	4" x 4" x .5" thick (10 cm x 10 cm x 1.2 cm)	5	A6209
HBHD6675	6" x 6" x .75" thick (15 cm x 15 cm x 1.9 cm)	5	A6210
<b>Dressings With Moisture-Retentive Film</b>			
HBRF4414	4" x 4" (10 cm x 10 cm)	10	A6209
HBRF2258	2.25" x 8" (5.7 cm x 20 cm)	10	A6210
HBRF4475	4" x 4.75" (10 cm x 12 cm) Island Dressing; <i>pad size 2" x 2.75" (5 cm x 6.5 cm)</i>	10	A6212
HBRF2650	2.5" (6 cm) Diameter Ostomy Dressing	10	A6209



Hydrate before use

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician (or properly licensed practitioner). Prior to the use of this device be sure to read the complete information on how to use this device, including Warnings, Cautions and Instructions for Use and all other package inserts supplied with the product.

For product questions, sampling needs, or detailed clinical questions concerning our products in the US call **860-337-7730**.

**Hydrofera**<sup>®</sup>

[www.hydrofera.com](http://www.hydrofera.com)

Hydrofera Blue Dressings  
Manufactured By:  
Hydrofera, LLC  
340 Progress Drive  
Manchester, CT 06042 USA

1. Instructions for Use
  2. Data on file 2013 – 130107-1CYT
  3. Data on file, EN 13726-1; 2002.
  4. Data on file 2017 - TR-00148
- \*HCPCS are for reference only and subject to change.

Hydrofera Blue and Hydrofera Blue READY are trademarks of Hydrofera, LLC.  
All other trademarks are the property of their respective owners.  
Made and printed in the USA.  
©2018 Hydrofera, LLC.  
002-0518