



Hydrofera Blue READY[®] and Hydrofera Blue TRANSFER[®] Dressings



Bioburden and Exudate Control at Your Fingertips

- Available with or without film backing
- Antibacterial properties to help manage bioburden
 - No hydration required prior to application
 - Up to 7-day wear time

When the Management of Bioburden Is a Priority



Broad-spectrum antibacterial protection¹

- Provided by two organic pigments, methylene blue and gentian violet
- Noncytotoxic²

Absorption

- Absorbs and retains bacteria-laden exudate

Available with or without waterproof film backing

Hydrofera Blue READY[®] and TRANSFER[®] Dressings Provide Broad-Spectrum Antibacterial Protection¹

Bacteria

- *Methicillin-resistant Staphylococcus aureus (MRSA)*
- *Staphylococcus aureus*
- *Staphylococcus epidermidis (coagulase negative)*
- *Pseudomonas aeruginosa*
- *Pseudomonas fluorescens*
- *Escherichia coli*
- *Enterococcus faecalis*
- *Vancomycin-resistant Enterococcus (VRE)*
- *Streptococcus pneumoniae*
- *Streptococcus pyogenes*
- *Klebsiella pneumoniae*
- *Proteus mirabilis*
- *Proteus vulgaris*
- *Enterobacter aerogenes*
- *Enterobacter cloacae*
- *Bacillus subtilis*
- *Yersinia enterocolitica*
- *Serratia marcescens*

Yeast

- *Candida albicans*
- *Candida krusei*
- *Candida glabrata*

Hydrofera Blue READY



- Waterproof film backing
- No secondary dressing required
- Exudate is absorbed and retained
- Up to 7-day wear time

Hydrofera Blue Transfer




- No film backing increases exudate management options
- Exudate is transferred into secondary absorbent dressing
- Suitable for use on partial and full thickness wounds
- Up to 7-day wear time

Performance Data

Absorption Comparison

Hydrofera Blue READY foam dressing has absorption comparable to other antimicrobial foam dressings, based on in vitro testing³




(Error bars represent one standard deviation)

Performance Data


In Vitro Fluid Transfer Study

Hydrofera Blue TRANSFER has a high rate of exudate transfer to a secondary dressing⁴




(Error bars represent one standard deviation)

Hydrofera Blue CLASSIC[®]



Product No	Product Size	Quantity/Box	HCPCS*
Standard Dressings			
HB2214	2" x 2" (5 cm x 5 cm)	10	A6209
HB4414	4" x 4" (10.2 cm x 10.2 cm)	10	A6209
HB4514	4" x 5" (10.2 cm x 12.7 cm)	10	A6210
HB6614	6" x 6" (15.2 cm x 15.2 cm)	10	A6210
HBT0906	9-mm Tunneling Dressing (1.2 g)	10 pair	A6215



Heavy Drainage Dressings			
HBHD4450	4" x 4" x 0.5" thick (10.2 cm x 10.2 cm x 1.2 cm)	5	A6209
HBHD6675	6" x 6" x 0.75" thick (15.2 cm x 15.2 cm x 1.9 cm)	5	A6210



Ostomy Dressing			
HBRF2650	2.5" (6.4 cm) Diameter Ostomy Dressing	10	A6209

Hydrofera Blue ComfortCel[®]




HBCC4550	4" x 5" x 0.5" (10.2 cm x 12.7 cm x 1.3 cm)	5	A6210
HBCC6675	6" x 6" x 0.75" (15.2 cm x 15.2 cm x 1.9 cm)	5	A6210

Hydrofera Blue READY[®]




HBRS2520	2.5" x 2.5" (6.4 cm x 6.4 cm)	10	A6209
HBRS4520	4" x 5" (10.2 cm x 12.7 cm)	10	A6210
HBRS8820	8" x 8" (20.3 cm x 20.3 cm)	10	A6211

Hydrofera Blue TRANSFER[™]



HBRT2525	2.5" x 2.5" (6.4 cm x 6.4 cm)	10	A6209
HBRT4050	4" x 5" (10.2 cm x 12.7 cm)	10	A6210
HBRT8080	8" x 8" (20.3 cm x 20.3 cm)	10	A6211

Hydrofera Blue READY-BORDER[®]



HBRB4040	4" x 4" (10.2 cm x 10.2 cm)	10	A6212
HBRB5959	5.9" x 5.9" (15 cm x 15 cm)	10	A6212
HBRB6080	6" x 8" (15 cm X 20.3 cm)	10	A6213

*REIMBURSEMENT DISCLAIMER: The Centers for Medicare & Medicaid Services (CMS) uses Palmetto GBA as the Pricing, Data Analysis, and Coding (PDAC) Contractor. The PDAC Contractor maintains the Product Classification List, located at https://www4.palmettogba.com/pdac_dmecs/. The Product Classification List is a searchable database containing products that have received HCPCS coding verification from PDAC and their corresponding HCPCS codes. Hydrofera LLC, however, makes no representations as to the accuracy of the information contained within the PDAC database, nor any representations as to whether its products are reimbursable under any government sponsored healthcare program and/or private-insurance program. It is the provider's sole responsibility, in consultation with the insurer, if necessary, to determine medical necessity, ensure coverage criteria is met, and submit appropriate codes, modifiers, and charges for services rendered.

1. Instructions for Use 2. Data on file 2013 - 130107-1CYT 3. Data on file, EN 13726-1; 2002. 4. Data on file 2017 - TR-00148
*HCPCS are for reference only and subject to change.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician (or properly licensed practitioner). Prior to the use of this device be sure to read the complete information on how to use this device, including Warnings, Cautions and Instructions for Use and all other package inserts supplied with the product.

For product questions and samples, please call
Customer Service at 1-860-337-7730.

Hydrofera Dressings manufactured by:
Hydrofera, LLC, 340 Progress Drive
Manchester, CT 06042, USA
www.hydrofera.com

Hydrofera, Hydrofera Blue, and the Hydrofera logo are registered trademarks of Hydrofera, LLC. ©2023
Hydrofera LLC is an Essity company.

14-001-001 Rev. 01 65142 R1 A24

MADE IN USA

Hydrofera[®]